

tactical care products

 Tactical medical equipment



Ladies and Gentlemen,
herewith you receive the new catalog

tactical care *products*

Our aim was to create a catalog specifically for the special requirements of medical equipment and equipment for armed forces and security agencies such as the police and customs, as well as the increased demands on the emergency services due to the current security situation.

In addition to products from well-known companies such as TacMed Solutions, Teleflex etc., we also offer an ever-expanding range of our own products under the 'Lifeguard Tactical Products' brand.

On the following pages you will find products on the following topics:

HAEMOSTASIS

SECURING THE AIRWAY AND VENTILATION

TREATMENT OF THORACIC INJURIES

HEAT PRESERVATION

BAGS, RUCKSACKS AND KITS

IMMOBILISATION

RESCUE AND TRANSPORT

AND MUCH MORE . . .

Benefit from our attractive prices, long opening hours, expert advice and our ability to deliver.

All products are available from well-stocked specialist retailers.

With kind regards
Your Lifeguard Team

We are certified
according to DIN EN ISO 13485:2021

Please address technical questions to:

Mr. André Damschen

Product Manager
Tel.: 0049 281-95283932
E-mail: andre.damschen@servoprax.de

Assortments and prices are subject to change. No liability for errors or printing errors.

Cover photo: Adobe Stock

© servoprax GmbH
Status: February 2025

10-point fixation system	13	Intraosseous infusion	5
		Intraosseous set	5
A		L	
AAJT-S Tourniquet	1	Laryngoscope blade	6
Abdominal Tourniquet	1	Laryngoscopes	6
Accessories for stretcher	14	Leg bag	9
Accessories scoop stretcher	13	Lifeguard Spineboard	12
Accessories spineboard	13	Lifeguard stretcher	14
Active thermal blanket	17	L-Mast anti-shock pants	16
AMP backpack	10	Loop belts	13
Angiocath	5	LTD stretcher rail	16
Anti-shock pants	16		
Aortic Tourniquet	1	M	
B		MLP leg bag	9
Backpack	10	Multi trauma bandage	2
Bag for SEL rescue sling	15		
Bandages	3	O	
Bandage scissors	11	Occlusis association	4
BD article	5	Oleas Bandage	2
Belt systems	13		
BLAST bandage	3	P	
C		Pelvic sling	13
Carrier	14		
Central access	5	Q	
Chassis for stretcher	14	Quick Breezer ventilation foil	6
Chest Seal	4, 18	Quicktrach	7
Clear Line II resuscitation bag	6		
Cold Blade Set	6	R	
Cold light laryngoscope set	6	Ready-Heat thermal blanket	17
Cold light spatula	6	Ready-Heat warming vest	17
Combat	2	Rescue blankets	17
Combat Tourniquet	1	Rescue board	12
Combi Carrier loop straps	13	Rescue scissors	11
Coniotomy set	7	Rescue sling	15
Cordura	9	RESQ Blade	6
Corps covering	19	Resuscitator bags	6
D		Robin Safety Boy scissors	11
Decompression needle	5		
Disposable Cold Blade	6	S	
E		ScalpelCric surgical coniotomy set	7
Emergency backpack	10	Scissors	11
Emergency bag	1, 8ff	Seldinger Coniotomieset	7
Emergency resuscitation bag	6	SEL rescue sling	15
Evacuation sheet	15	Sof TT-W	1
Evacuation stretcher	15	Spare filling emergency bags	8ff
Exercise leg	18	Spatula for laryngoscopes	6
Eye patch	11	Special cannula	5, 7
Eye Shield	11	Special scissors	11
F		Spider strap	13
FAP bag	8	Spinal board	12
Fast Responder	5	Spineboard	12
Filling assortment CFR	10	Splint material	16
Filling assortment MLP leg bag	9	Stainless steel laryngoscope blade	6
First Aid backpack	10	Stretcher	14f
First Aid bag	1, 8ff	Stretcher chassis	14
Foxtrot Litter Stretcher	15	Stretcher rail	16
Frame for stretcher	14	Surgical coniotomy	7
G		Surgicric	7
Gauze bandages	2, 18	T	
Gluing tape	2	TacShears	11
H		Tape	2
Heat Vest	17	TCP Bag	1
Hemostasis	1ff, 18	Thermal blanket	17
Hypothermia protection	17	Thoracic drainage	5
I		Thoracic occlusion plaster	4
IFAK	9	Tourniquet	1
IFAK filling assortment	8	T-POD Combat pelvic sling	13
Infusion set	5	Trainingsbandage	18
Intraosseous cannula	5	Trauma bandage	2
		Trauma training	18
		W	
		Wound Cover	2f
		Wound Packing Trainer	18



SOF Tourniquet

The newest evolution of the SOF Tourniquet sets the benchmark for prehospital tourniquets. Its innovative design incorporates end-user feedback from specific experiences and features most important to them, creating a tourniquet that is stronger, faster, and smoother than ever before.

Updates

- » Easy-pull action for smooth, quick self-application
- » Custom engineered materials for increased strength and durability
- » Lighter weight for improved portability

Features:

- » Completely redesigned („1.5“ width) webbing maximizes compression during the application process. The new performance material glides through the Rugged Buckle for quick, snag-free operation and minimizes loss of pressure over time after application
- » New rugged buckle
- » Slack Indicator Wedge: stitched into the Performance Compression Band below the windlass, the contrasting wedge offers a visual confirmation you have pulled all excess slack from the band.
- » Reinforced Tourniquet Retention Assistance Clip

Colour	SU	REF	Price SU
black	1 piece	N7 1007	39,50
orange	1 piece	N7 1007-O	39,50
blue	1 piece	N7 1007-B	39,50
tan	1 piece	N7 1007-T	39,50
red	1 piece	N7 1007-R	39,50



TCP-Pouch

The TCP-Pouch (Tourniquet Carrier Pouch) is a small pouch to securely carry your Tourniquet. It's MOLLE-compatible and therefore can be attached to your vest or any other gear with MOLLE. Empty.

Technical data: Dimensions: 20 x 6 x 4 cm, Weight: 50 g

Colour	SU	REF	Price SU
black	1 piece	N7 TCP01	9,95
olive-green	1 piece	N7 TCP03	9,95



AAJT-S Tourniquet

Non-compressible hemorrhages are the cause of death for tens of thousands annually. Many difficult injuries in the abdomen, axilla, pelvis and inguinal areas require immediate care by first responders, paramedics, EMS professionals, emergency physicians, police officers and firefighters to help stabilize accident victims. The ability to control difficult to manage hemorrhages quickly and easily is imperative for patient survival until they can reach definitive care.

The Abdominal Aortic and Junctional Tourniquet-Stabilized (AAJT-S) device by Compression Works, was created to provide a solution to control difficult to manage junctional hemorrhage and pelvic bleeding. The AAJT-S is a multi-functional, externally applied truncal and junctional tourniquet that utilizes a wedge-shaped bladder to press into the body to occlude blood flow through all vasculature under the site of application.

The AAJT-S is cleared by the FDA for use to control difficult bleeding situations in the abdomen, axilla, inguinal area and for pelvic fractures. The AAJT-S device has been used successfully in military combat and in civilian pre-hospital settings and is the only junctional tourniquet proven to save human life in upper and lower junctional bleeding. The device also provides pelvic fracture stabilization and has been shown to be equivalent to Zone 3 REBOA.

NSN: 6515-01-616-4999

Technical data: Dimensions: 23 cm x 20 cm x 5 cm, Weight: 425 g

Type	SU	REF	Price SU
AAJT-S Tourniquet	1 piece	N7 1300	499,00
AAJT-S Tourniquet Demo-Unit	1 piece	N7 1300T	299,00



Lifeguard Multi Trauma Bandage

The Multi Trauma Bandage has been developed to provide a large number of different injury patterns with a bandage. The pressure point is made for a better application of a pressure bandage, the twofold wound dressing can be used for extensive injuries (e.g. in abdomen) or else for simultaneous treatment of entry and exit wounds.

The use of the Multi Trauma Bandage is intuitive, as it is applied like each normal dressing.

Velcro strips inside the bandage prevent slipping, loosening and unintentional rolling off.

The Multi Trauma Bandage is single sterile vacuum-packed.

Dimensions	OC	SU	REF	Price SU
10 cm x 4,5 m (elastic)	100 pcs.	1 piece	N7 MTB4	7,10
15 cm x 4,5 m (elastic)	100 pcs.	1 piece	N7 MTB6	7,70

OLAES Bandage

No other trauma bandage of this size is as versatile as the OLAES Bandage. The bandage can be used to dress minor wounds, for pressure dressings, occlusive dressings or treat severe wounds with the occlusive plastic sheet and dressing pad. Equally useful for a qualified paramedic or untrained first-aid helper, the OLAES Bandage should be in all medical and first-aid kits.

Every first-aid-er could use the OLAES Bandage, as it may be applied like any normal elastic bandage. There are no hooks or clips to complicate usage.

Pressure can be applied to the wound with the help of the pressure bar and thus control haemorrhaging. Velcro strips hold the OLAES Bandage in place during application. Each bandage is sterile and singly vacuum-packed.

NSN: 6510-01-558-3342 - REF: N7 1001

NSN: 6510-01-631-7893 - REF: N7 1001-F

NSN: 6510-01-558-3443 - REF: N7 1002



Type	Dimensions	Colour	SU	REF	Price SU
rolled	10 cm x 3 m	beige	1 piece	N7 1001	12,45
sheets	10 cm x 3 m	beige	1 piece	N7 1001-F	12,45
rolled	15 cm x 3 m	beige	1 piece	N7 1002	13,25



Lifeguard Compact gauze bandages sterile

White, with Z-fold, 6-layers, 10 cm x 3.75 m.

Sterile, vacuum packed.

Type	SU	REF	Price SU
sterile	1 piece	J8 10005	2,10



Lifeguard Combat Tape

The Combat Tape is a universally applicable, strong adhesive tape that gives first responders additional security and support when gauze and elastic bandages are not enough. The roll is small and handy and therefore fits in every IFAK, but with the size of 5 cm x 2.5 m it offers enough tape for a wide range of applications.

Areas of application:

- » Emergency services
- » Emergency medical services
- » Military
- » Police
- » Disaster relief

Technical data: Version: Color silver, Dimensions open: 5 cm x 2.5 m, Dimensions packed: 5 x 3.9 cm, Weight: 40 g

Dimensions	SU	REF	Price SU
5 cm x 2,5 m	1 piece	N7 CT1	2,65





Lifeguard Abdominal Trauma Bandage

The Abdominal Trauma Bandage ATB was designed for the treatment of extensive injuries to the abdomen as well as for traumatic amputations and burns. For this purpose, the bandage has a non-adhesive wound dressing of 50 x 50 cm, which can be placed on or around the injured region and securely fixed with the elastic bandage.

Velcro strips at intervals on the bandage ensure a secure hold and also that the bandage cannot unroll unintentionally.

The ATB is individually sterile and space-saving vacuum packed.

Dimensions	OC	SU	REF	Price SU
15 cm x 4,5 m (elastic)	100 pcs.	1 piece	N7 ATB6	10,00



BLAST Bandage

The BLAST Bandage was developed to enable traumatic amputations, burns and large wounds to be dressed quickly and with little fuss. The wrapped dressing is 10cm wide, the unwrapped dressing can cover a wound of approximately 50x50cm. The non-adherent wound pad can be to open wounds and limbs and held in place using the elastic bandage. The BLAST Bandage is large enough to cover the entire back or chest of large casualties.

The BLAST Bandage offers a removable occlusive layer 47.5 x 47.5 cm which can be used to cover abdominal wounds and thus prevent loss of heat and moisture. Each bandage is sterile and singly vacuum-packed.

NSN: 6510-01-586-4314

Dimensions	Colour	SU	REF	Price SU
50 x 50 cm	white-beige	1 piece	N7 1000	14,00



Lifeguard Easy Bandage

The Easy Bandage is used for quick wound care. Velcro strips inside the bandage prevent slipping, loosening and unintentional rolling off, whereby a contamination of the bandage is prevented as well.

The Easy Bandage is available in 2 different widths.

Dimensions	OC	SU	REF	Price SU
10 cm x 4,5 m (elastic)	240 pcs.	1 piece	N7 EB4	2,65
15 cm x 4,5 m (elastic)	240 pcs.	1 piece	N7 EB6	2,99



Control Wrap Bandage

Using the Control Wrap bandage, wounds may be dressed quickly and efficiently. The Velcro strip ensures that the bandage does not slip or unroll whilst being applied, thus reducing the risk of contamination or loss of pressure on the wound and also saves time.

Dimensions	Colour	SU	REF	Price SU
10 cm x 1,80 m	beige	1 piece	N7 1003	4,65
15 cm x 1,80 m	beige	1 piece	N7 1004	4,95



HALO SEAL

First aid for open chest wounds

HALO SEAL is one of the original chest seals, intended to treat penetrating chest injuries. HALO SEAL acts as a mechanical barrier, adheres in the presence of moisture, dirt, blood or other residues.

Product Features

- » Durable packaging
- » Large pull tab for easy removal
- » No vent to occlude
- » No latex used in the manufacturing of the product

Product Specifications

- » 6-year shelf life
- » Two seals, treat entrance and exit wounds
- » Performs in temperature range of -34° – 60° C.
- » Dimensions: ca. 16.5cm x 14cm
- » Manufactured in the USA

Pack with 2 pcs.

NSN: 6510-01-581-0553

Type	Dimensions	SU	REF	Price SU
double package	14 x 21 cm	1 pack	N7 1163	23,00



Rhino Rescue Chest Seal

The Rhino Rescue Chest Seal is an advanced occlusive dressing designed to treat open chest wounds. Penetrating injuries to the chest that result in open chest wounds may lead to tension pneumothorax, the second most common cause of preventable death on the battlefield.

The Rhino Rescue Chest Seals are specifically designed with a hydrogel adhesive strong enough to adhere to skin – regardless of sweat, blood, or hair.

Features:

- » Used by Military, Police, Fire, EMS and First Responders
- » with valves
- » Large Release Tab for Easy Application and Removal
- » Effective in Extreme Cold and Heat
- » Transparent Backing for Easy Placement Over the Wound
- » Strong Hydrogel Adhesive
- » Conforms and Adheres Easily to All Body Surfaces
- » Provides total occlusion – even with excessive blood, dirt, sand, hair, or heavy perspiration
- » Configured for ease of use in low-light environments
- » Latex-free
- » Durable, puncture-resistant packaging
- » Integrity and function maintained when stored or carried in folded configuration

Dimensions	SU	REF	Price SU
15 x 9,5 cm	1 piece	N7 1162	15,50



FAST1 Intraosseous Infusion System

The FAST1 Intraosseous Infusion System is designed to establish a sternal intraosseous access route to patients when standard venous access is not possible. Using the sternal IO access route, providers can rapidly deliver lifesaving fluids and medications to facilitate emergency resuscitation through the bone marrow space. FAST1 can be effectively used during CPR and when a cervical spine collar is in place.

- » Sternum insertion site
- » Flow rates of up to 80mL/min (gravity drip)
- » Depth control (no need to choose needle length)
- » Low Profile Tubing (infusion tubing designed to be secure and not dislodge)
- » Small footprint 21.59 x 20.32 x 3.30 cm
- » Lightweight 127.6 g
- » Manual insertion – no batteries required
- » Indicated for patients aged 12 or older

Type	SU	REF	Price SU
FAST1 Intraosseous Infusion System	10 pcs.	I4 8001	3.360,00



Complete thoracic drainage tray

Set completely for applying the drainage in pragmatic manner, in an ambulance or rescue van within emergency situations.

Contents:

- 1 XRO trocar drain
- 1 fenestrated drape 60 x 50 cm
- 2 rectangular cups for detergent and disinfecting agents
- 4 pads
- 8 gauze dressings 10 x 10 cm
- 2 dressing forceps
- 2 syringes 10 ml
- 1 needle for drug aspiration 38 x 1.2 mm (18G)
- 1 needle for skinwheel 15 x 0.5 mm (25G)
- 1 needle for local anaesthesia 38 x 0.8 mm (21G)
- 1 needle for exploratory puncture 90 x 1.2 mm (18G)
- 2 suture needles
- 1 scalpel
- 1 dual suction valve
- 1 plastic bag with connector for collecting the drained liquid
- 1 coil of adhesive bandage 75 mm wide
- 1 pair of gloves

Type	SU	REF	Price SU
Complete tray for thoracic drainage	1 set	N2 639-28	116,50



Angiocath Special Cannula 14 G

Special cannula, 83 mm long, for relief of a tension pneumothorax. Furthermore, can be used with obese and muscular patients. The indwelling catheter is much harder and more stable than a conventional IV catheter.

Size	OC	SU	REF	Price SU
14 G x 3 1/4	50 pcs.	10 pcs.	I4 382268	48,50
14 G x 5 1/4	50 pcs.	5 x 10 pcs.	I4 382269	0,00



T-Pack Special cannula 14 G

Special cannula for relief of a tension pneumothorax. This cannula is 83 mm long and thus sufficient for obese and muscular patients. The indwelling catheter is much harder and more stable than a conventional IV catheter.

Type	Size	SU	REF	Price SU
T-Pack cannula	14 G x 3 1/4	1 piece	I4 TM303	20,95



RESQ-BREEZER RESUSCITATOR BAG

Clearline II Resuscitator

Single Use-Resuscitator-Sets made of PVC, it includes a mask, an O₂-Reservoir bag and an O₂-Connection tube.

- » Latex free
- » Single packed in a poly bag
- » With all-in-one intake valve

Clear Line II	Adult	Child	Infant
	E5 04010-S	E5 04020-S	E5 04030-S
Material Corpus	PVC	PVC	PVC
Volume	1500 ml ± 200 ml	550 ml ± 200 ml	280 ml ± 100 ml
Disinfectable/ autoclavable	No	No	No
Pressure valve	40 - 60 cm H ₂ O	35 - 50 cm H ₂ O	30 - 45 cm H ₂ O

Type	Size	Volume	SU	REF	Price SU
Adult	Size 5	1500 ml ± 200 ml	1 Set	E5 04010-S	14,70
Child	Size 3	550 ml ± 200 ml	1 Set	E5 04020-S	14,70
Infant	Size 0	280 ml ± 100 ml	1 Set	E5 04030-S	14,70



RESQ-BLADE

RESQ-Blade Cold Blade ECO-LED-SET

Comprises of: each of 1 disposable metal blade: type Macintosh in sizes 2, 3, and 4 and type Miller Gr. 0 and 1, with disposable cold-light LED-battery handle in a sturdy, blue roll pouch.

Type	SU	REF	Price SU
Cold Blade ECO-LED-Set	1 set	I5 091LED	29,50



QUICK-BREEZER

Quick-Breezer ECO Face shield

Rapid resuscitation

The simple resuscitation face shield without a mouthpiece for protecting helpers and patients. Comes in an alu blister and provides:

- » Safety and hygiene during mouth to mouth and mouth to nose resuscitation
- » Very easy handling with the help of a brief instruction
- » With minimum space requirements

Type	SU	REF	Price SU
Quick Breezer ECO face shield, packed singly in alu-blister; without mouthpiece	1 piece	N4 4000	1,60



Quick Trach tracheotomy set

The tracheotomy set with an anatomically shaped cannula.

It creates an air passage simply and quickly. The plastic cannula with an inner diameter of 4 mm for adults and 2 mm for children enables sufficient ventilation using a resuscitator bag. It can be used with the mouth or with a conventional respiration device. The conical needle point can be inserted without prior incision – fast, small stoma, reduced risk of bleeding and swabbing. The stopper prevents excessively deep penetration and the easy-glide tip enables clear identification of the trachea. The 15 mm connector provides quick attachment to a resuscitator bag or respiration device, with or without a connecting tube. The holding plate guarantees secure fixation using the foam neck strap. Sterile and complete with foam neckband, syringe and connection tube.

Type	Inner diameter	SU	REF	Price SU
Child, in blister	2 mm	1 piece	N2 3004002	102,80
Child, plastic tube, in blister	2 mm	1 piece	N2 3004902	102,80
Adult, in blister	4 mm	1 piece	N2 3004004	102,80
Adult, plastic tube, in blister	4 mm	1 piece	N2 3004904	102,80



ScalpelCric surgical set for cricothyrotomy

ScalpelCric is a surgical set for cricothyrotomy according to the scalpel technique. The set has been developed to match the latest recommendations from the DAS (Difficult Airway Society) 2015 guidelines. Based on the simple description of „stab, twist, bougie, tube“ the ScalpelCric is meant to facilitate the execution of the scalpel technique. Each component is designed to enhance the safety and efficacy of the technique.

The set is packed in a specific pouch and the individual components are organised in the order of use.

Content: Scalpel fig.10, 14 Fr bougie (40 cm), 6.0 mm I.D. cuffed tube, Extension tubing, Syringe 10 ml, Necktape

Type	SU	REF	Price SU
ScalpelCric surgical set for cricothyrotomy	1 set	N2 30084471	92,82



Surgicric III Surgical Cricothyrotomy Set of the Seldinger technique.

The set Surgicric III contains a unique introducer needle for the optimal skin puncture with temporary oxygenation via the 15 mm connector and the implementation of the Seldinger technique.

Type	SU	REF	Price SU
Surgicric III Surgical Cricothyrotomy Set	1 set	N2 30082271	124,85



Lifeguard FAP-Pouch

The FAP-Pouch is used for an individual first aid kit for penetrating wounds. This compact pouch is MOLLE compatible and can be worn on the assault vest/body armor attachment. It offers a small Velcro pad for identification. Empty.

The FAP is available in 2 colours.

Technical data: Size: 21 x 13.5 x 4.5 cm, Weight: 160 g

Colour	SU	REF	Price SU
black	1 piece	N7 FAP01	26,30
coyote (beige brown)	1 piece	N7 FAP02	26,30



Basic Filling for IFAK (Consistent with FAP-Bag)

Quantity	Type	REF
1 pc	SOFTT-W tourniquet	N7 1007
1 pc	Multi Trauma Bandage 4, 10 cm x 4.5 m	N7 MTB4
1 pc	Quick-Breezer ECO ventilation barrier	N4 4000
1 pc	Emergency blanket green-silver	H1 16904
1 pc	NasoSafe tube CH28	I7 038-95-065
1 pc	TacShears Rescue shears	N7 002
1 pair	SoftHand nitrile-gloves L	H9 7002
1 pc	M-Silk plaster, 2.50 cm x 5 m	J8 702452

Type	SU	REF	Price SU
Filling Base for IFAK	1 piece	N7 BASIC	59,95



Expanded IFAK Filling (Consistent with FAP-Bag)

Quantity	Type	REF
1 pc	SOF tourniquet	N7 1007
1 pc	Multi Trauma Bandage 4, 10 cm x 4.5 m	N7 MTB4
1 pc	HemCon ChitoGauze XR PRO 7.5 cm x 3.7 m	N7 1090
1 pc	Compressed gauze, sterile	J8 1000
1 pc	Quick-Breezer ECO ventilation barrier	N4 4000
1 pc	Emergency blanket green-silver	H1 16904
1 pc	NasoSafe tube CH28	I7 038-95-065
1 pc	Optilube lube 5 gr	I7 4343
1 pc	TacShears Rescue shears	N7 002
1 pair	SoftHand nitrile-gloves L	H9 7002
1 pc	M-Silk plaster, 2.50 cm x 5m	J8 702452

Type	SU	REF	Price SU
Expanded IFAK Filling Range	1 piece	N7 ADVANCE	129,00





Lifeguard MLP-Leg Rig Puch

The MLP is a robust Leg Rig Pouch, that offers enough space for an extensive content of first aid equipment for several kinds of wounds. A rescue scissor can be placed on the backside of the pouch and on the front side is a Velcro pad for identification. Empty.

Technical data: Size: 22 x 20 x 8 cm, Weight: 350 g

Colour	SU	REF	Price SU
black	1 piece	N7 MLP01	35,90
coyote (beige brown)	1 piece	N7 MLP02	35,90



Lifeguard IFAK Pouch with Laser-Cut MOLLE, Cordura

The Advanced IFAK Pouch is a newly designed small emergency bag to store your individual first aid equipment. Made from robust Cordura (500D), it meets even the toughest requirements for such a bag and with the well thought-out interior division you can fill the bag clearly.

With the help of the laser-cut MOLLE on the front, additional small pockets such as our tourniquet carrier can be attached. If you do not use the item, you have a largely smooth surface that prevents you from getting caught in the terrain.

Technical data: Material: Cordura 500D,

Open dimensions: 20 cm x 16 cm x 5 cm, Laser-cut MOLLE, Weight: 160 g

Colour	SU	REF	Price SU
black	1 piece	N7 AIP01	30,50



Filling according to First Responder

(suitable for MLP-leg bag)

Quantity	Type	REF
1 pc	SOF tourniquet	N7 1007
1 pc	MTB Multi Trauma Bandage 4, 10 cm x 4.5 m	N7 MTB4
1 pc	Blast bandage 50 x 50 cm	N7 1000
1 pc	HemCon ChitoGauze XR PRO 7.5 cm x 3.7 m	N7 1090
1 pc	Quick-Breezer ECO ventilation barrier	N4 4000
1 pc	triangular bandage, cotton	H1 15804
1 pc	Emergency blanket green-silver	H1 16904
1 pc	TacShears rescue shear	N7 002
1 pair	Softhand nitrile-gloves L	H9 7002
1 pc	M-Silk plaster, 2.50 cm x 5 m	J8 702452

Type	SU	REF	Price SU
Filling following the equipping of the first-aid Alpha	1 piece	N7 EEH-A	139,00



LIFEGUARD
TACTICAL PRODUCTS



Lifeguard AMP Assault Medic Pack

The AMP (Assault Medic Pack) has been designed for short time missions and offers enough space for all accessories needed for the rescuer. Inside the backpack offers a panel with elastic loops and straps for fixing small items (Decompression needles for example) and 6 mesh pockets (most of them removable) for a good storage and overview of your equipment. With help of the two loops on the top the backpack can be hung and used as a panel. Furthermore, the AMP can be added to the MMP-backpack by zipping both together, giving you a complete solution for transportation of your expanded medical equipment. Empty.

Technical data: Dimensions: 46 x 30,5 x 10 cm, Weight: 1.2 kg

Colour	SU	REF	Price SU
black	1 piece	N7 AMP01	120,80
coyote (beige brown)	1 piece	N7 AMP02	120,80



Filling Expanded for CFR (Consistent with AMP-Backpack)

Quantity	Type	REF
1 pc	SOF Tourniquet GEN 4	N7 1007
3 pcs	MTB Multi Trauma Bandage 4, 10 cm x 4.5 m	N7 MTB4
1 pc	MTB Multi Trauma Bandage 6, 15 cm x 4.5 m	N7 MTB6
3 pcs	Easy bandage 10 cm	N7 EB4
1 pc	ATB Abdominal Trauma Bandage 6, 15 cm x 4.5 m	N7 ATB6
1 pc	HemCon ChitoGauze XR PRO 7.5 cm x 3.7 m	N7 1090
2 pcs	Compressed gauze	J8 1000S
1 pc	NasoSafe tube CH28	I7 038-95-065
1 pc	Optilube lube 5 gr	I7 4343
1 pc	ScalpelCric Surgical Coniotomy Set	N2 30084471
2 pcs	Quick-Breezer ECO ventilation barrier	N4 4000
1 pack	HALO Chest Seal 14 x 21 cm, double packing	N7 1163
2 pcs	T-Pack decompression needle	I4 TM303
2 pcs	E-Bone Splint	N5 0250
1 pc	TacShears rescue shear	N7 002
1 pair	Softband nitrile-gloves L	H9 7002

Type	SU	REF	Price SU
Filling expanded for CFR	1 piece	N7 CFR	499,00



Lifeguard Eye shield, eye protection made of polycarbonate

The LIFEGUARD TACTICAL PRODUCTS Eye Shield is a convex rigid eye shield designed to protect the injured eye from further external influences that could cause further damage to the eye.

It complies with current military and civilian guidelines for the treatment of traumatic eye injuries in the field.

Characteristics:

- » Made of rigid, shatterproof polycarbonate
- » Meets CoTCCC treatment recommendation for ocular trauma
- » Air circulation possible through holes in the eye shield
- » Suitable for every eye
- » Translucent orange color

Areas of application:

- » Emergency services
- » Emergency medical services
- » Military
- » Police
- » Disaster relief

NSN: 6515-01-647-5456

Technical data: Open dimensions: 7.5 x 6.2 x 0.3 cm,
Dimensions packed: 7.5 x 6.2 x 0.3 cm, Weight: 7.5 g

Type	Colour	SU	REF	Price SU
Eye Shield, eye protection made of polycarbonate	orange	1 piece	N7 ES1	2,85



TacShears, Rescue Shears, Black

Rescue shears, which has been designed for tactical purposes, with black titanium coated blades for a particularly good cut, with plastic handles and blunts ends.

Type	Colour	Dimensions	SU	REF	Price SU
TacShears Rescue shears	black	18 cm	1 piece	N7 002	7,92



Robin Safety-Boy emergency scissors

These scissors are multifunctional:

Ideal for cutting through safety belts and clothing and the sturdy steel makes it possible to break through glass windows. Without spring.


Type	Colour	Dimensions	SU	REF	Price SU
Robin Safety-Boy Emergency scissor	chrome	23 cm	1 piece	N9 4900	71,00



Forma emergency scissors

These multifunctional scissors are ideal for cutting through belts and clothing and the sturdy Solingen steel makes it possible to break through glass such as car windows. The massive scissors come in a special holder which can be attached wherever necessary so that the scissors are always at hand.

Type	Dimensions	SU	REF	Price SU
Forma emergency scissors	23 cm	1 piece	N9 2010	46,50



Lifeguard Spineboards

Made of sturdy polyethelene, without joins or seams, therefore watertight.

16 large handgrips around the perimeter to make it easier to transport the patient. The surface is moulded for comfort, the foot end is slightly angled to simplify the rescue of a patient. The Spineboard conforms with EN 1865-1 and is suitable for water rescue, as well as being 100 % radioluscent.

Available in the colours: Olive green (tactical green) and black (tactical black)

Suitable for Speedclip harness attachment due to 4 continuous Speedclip bars.

Both Boards are only available as the Speedclip types.

Technical data: Length: 1,860 mm, Width: 406 mm, Height: 55 mm, Weight: 8.8 kg,
Max. load: 200 kg, Buoyancy during water rescue: up to 135 kg load

Type	SU	REF	Price SU
Lifeguard Spineboard, green (tactical-green)	1 piece	N5 0300-TG	259,00
Lifeguard Spineboard, black (tactical-black)	1 piece	N5 0300-B	259,00



Lifeguard Spineboard spider strap

10-point fixation system made of black Cordura-Nylon, with a central strap and 4 transverse straps, 3 of which are adjustable to enable quick application and optimise fixation. The foot and shoulder straps are fixed. Includes storage and transport bag.

Adjustable, according PHTLS-directive.

Type	Colour	SU	REF	Price SU
10-point fixation system, including bag	black	1 piece	N5 ST02035-B	69,00



Lifeguard Head Fix head fixation

Universal head fixation for shovel stretcher, spineboard and Combi Carrier, completely washable made of PVC-free material.

Soft foam with Velcro straps for fixation. The ears remain uncovered, X-ray compatible.

Colour	SU	REF	Price SU
blue	1 piece	N6 SH00111C	79,00
black	1 piece	N6 SH00121B	79,00



T-POD Combat Pelvic stabilisation device black

The T-POD Combat-Pelvic stabilisation device is a rapid and effective device for treating pelvic trauma, pain, and bleed-out while preventing loss of life. With its patented circumferential tightening mechanism, T-POD Responder is the only external fixation device shown in independent clinical studies to be equivalent to a surgical procedure in achieving stabilisation, while also mitigating patient discomfort. T-POD is used extensively by civilian first responders in motor vehicle accidents and by medics on the battlefield, especially in incidents involving improvised explosive devices.

T-POD Advantages

- » Complete: specially engineered pulley system ensures even pelvic compression and stabilisation
- » Easy: single-handed tightening with no straining
- » Adjustable: one size fits most physiologies, easily trimmed for custom fit; two can be secured together for obese patients
- » Flexible: non-metal fabric belt
- » Compatible: 100% radiolucent, X-ray, CT, and MRI compatible
- » Safe: latex-free
- » Non sterile - non reusable

Type	Colour	SU	REF	Price SU
T-POD Combat Pelvic stabilisation device	black	1 pcs.	N7 1100-C	159,00



Lifeguard Griffin stretcher

The „Griffin“ is a tactical four-foldable stretcher, made of 38 mm thick anodized aluminum tubes for braving the most adverse circumstances on the battlefield. The black cover is made of synthetic net-material which is easy to scrub. According to its purposes the stretcher is light and compact with a total weight of 9 kg and because of that absolutely rugged. With the transport backpack (included in the standard delivery) the stretchers are easy to stow and carry to the operation area.

Technical data: Technical details, Total weight: 9 kg, Size (stretching out): 2290 x 570 x 150 mm, Size (folded): 530 x 180 x 180 mm, Max. load: 220 kg

Type	SU	REF	Price SU
Griffin stretcher	1 piece	N7 G2000	573,30



Carrier EL3000

The chassis is a real all-rounder. With its extra large tires, it has been developed specially, to transport injured persons or material easily and quickly through impassable terrain. Additionally, the unique design allows, that an injured person can be transported by two or even one person, instead of 4 persons as before. This makes sense everywhere, where many injured must be supplied at the same time (Disaster areas, trouble spots etc.) or the rescuers must be protected themselves.

Another advantage of the carrier is, that it can also be used in an emergency as an additional camp bed or examination couch, by simply keeping the patient on it.

Thanks to its compact folding design, it can be stowed away in a space-saving manner when not in use. It is compatible with the most common stretchers as well as spineboards and should not be missing anywhere, where rough terrain or long distances have to be overcome (medical services, water rescue, mountain rescue, airfields etc.).

NSN: 6530-15-192-9999 - olive green

Technical data: Length closed/open: 52 cm / 158 cm, Width: 62 cm, Height closed/open: 84 cm / 86 cm, Weight: frame 20 kg, tyre 4 kg, Max. load capacity: 200 kg

Type	Colour	SU	REF	Price SU
Carrier EL3000	olive-green	1 piece	N7 EL3000	1.440,00



Foxtrot Litter DA

The Foxtrot Litter is an effective solution for rescuing injured persons, is simple and quick to use and small enough for one person to carry. Useful regardless of whether the situation requires a quick response or is of a longer duration. The compact size allows the Foxtrot Litter to be added to a backpack and therefore always at hand should it be required. The litter is rolled out and the injured person placed onto it and quickly transported without any delay caused by complicated fastening or fixation. The integrated ropes enable the Foxtrot Litter to be pulled along.

Features:

- » Light and compact
- » Enables transport of injured through narrow passages
- » The injured can be stood up to negotiate tight corners in stairwells
- » Draglines quickly unpacked
- » Several carry handles
- » Safety belts quick to apply
- » May be dragged or carried
- » The injured may be treated whilst being transported

NSN: 6530-01-571-5667

Technical data: Dimensions: 212 cm x 46 cm,
Size when packed: 53 cm x 16 cm, Weight: 1990 g, Max. load: 200 kg

Type	SU	REF	Price SU
Foxtrot Litter DA	1 piece	N7 LG3301	285,00



SEL-Litter (Small Evacuation Litter)

THE SEL-Litter is a lightweight and compact litter for evacuation of casualties from the danger zone. It's equipped with a bag for the feed as well as a fixation strap, so that the casualty can easily and fast be transported. The SEL has 8 grab points.

Because of the also available pouch the SEL can be stored easily via MOLLE on your combat vest or any other equipment.

Technical data: Weight 400 g, Max. load: 159 kg

Type	Colour	Dimensions	SU	REF	Price SU
SEL-Litter	black	180 x 70 cm	1 piece	N7 7520	49,40



Carrying Pouch for SEL-Litter

Small Carrying Pouch for the SEL Litter. Can be adopted to your equipment via MOLLE.

Technical data: Size: 30 x 22 cm

Colour	SU	REF	Price SU
black	1 piece	N7 7521	18,40
coyote (beige brown)	1 piece	N7 7522	18,40
olive-green	1 piece	N7 7523	18,40



LIFEGUARD®
EMERGENCY PRODUCTS



Lifeguard Anti-shock trousers

The revolutionary L-MAST provides safe, effective pressures without the use of inflatable bladders. Made from high-quality durable neoprene. L-MAST is simple to use and easy to clean. The non-inflatable L-MAST maintains constant pressure at varying altitudes and temperatures making it ideal for air rescue operations. L-MAST is lightweight and compact.

Colour: black

NSN: 6515-17-112-7790

Technical data: Dimensions: 38 x 28 x 15 cm (folded), Weight: 1.5 kg

Type

Lifeguard Anti-Shock Trousers, bag included

SU REF Price SU

1 piece N5 31801 203,49



LIFEGUARD®
EMERGENCY PRODUCTS

Lifeguard LTD Traction Device

The LTD is a high quality traction device which alleviates pain caused by fracture whilst. It prevents the fracture from damaging neighbouring arteries or nerves. Includes instructions and a bag.

- » Lightweight and space-saving
- » Can be used under and over MAST
- » The leg must not be raised to apply the LTD
- » Splinting a fractured humerus is possible

Type

LTD Traction Device

SU REF Price SU

1 piece N7 3000BK 198,50



Ready-Heat II Thermal blanket

This thermal blanket protects against hypothermia by actively creating heat. Made of a non-woven fabric interwoven with 12 chemical elements which reacts with the air to create warmth if the blanket is removed from the vacuum packaging. The blanket heats to ca. 38 °C within 10 minutes (possible max. temperature 52 °C) and retains this temperature for ca. 8 hours.

NSN: 6532-01-538-0759

Technical data: Dimensions: 32.25 x 40 cm (packed), 86 x 122 cm (open), Weight: 0.68 kg

Type	SU	REF	Price SU
Ready-Heat II Thermal blanket	1 piece	N6 RH2	38,50



Ready-Heat Thermal Blanket

Blanket with heat pads that warms to a minimum of 40 °C within 15 - 20 minutes of breaking the package seal. The blanket will maintain that temperature for 10 hours. Used in the prevention of hyperthermia and temperature management.

Suitable as an evacuation litter and able to carry a person who weighs up to 113 kg.

NSN: 6532-01-525-4062 - 86 x 122 cm

NSN: 6532-01-525-4063 - 86 x 152 cm

Type	Dimensions	SU	REF	Price SU
Ready Heat Air activated blanket, with 4 heat pads	86 x 122 cm	1 piece	N6 MD1	54,60
Ready Heat Air activated blanket, with 6 heat pads	86 x 152 cm	1 piece	N6 LG1	57,50



Ready-Heat Air activated blanket small

Air activated blanket that warms to a minimum of 40 °C within 15 minutes of breaking the package seal. The blanket will maintain that temperature for 10 hours.

Has two adjustable hook and loop straps which allow attachment to injured extremities. Can also be used as an IV fluid warmer.

NSN: 6532-01-538-1525

Technical data: Measurement: 41 x 46 cm

Type	SU	REF	Price SU
Ready Heat Air activated blanket	1 piece	N6 SM1	20,95



Air activated vest

Air activated vest that warms to a minimum of 38 °C within 10 - 15 minutes of breaking the package seal. The vest will maintain that temperature for 8 hours.

Overhead one size fits all style with two adjustable hook and loop straps.

Technical data: Universal size.

Type	SU	REF	Price SU
Ready-Heat Air activated vest	1 piece	N6 V01	30,50



Simu-Seal

The Simu-Seal has been designed to promote the use of medical chest seals and familiarise Military and Law Enforcement personnel with the indications and use of medical chest seals in the setting of training. This device is both INEXPENSIVE and EASY TO REMOVE from even densely hair covered skin with minimal or no discomfort. Made in the U.S.A. from high quality, latex free, medical grade materials.

Type	SU	REF	Price SU
single package	1 pack	N7 1167	5,65
double package	1 pack	N7 1168	8,70



HemCon ChitoGauze Trainer

The HemCon ChitoGauze Trainer is used for training in using the Hem-Con ChitoGauze bandages.

Type	SU	REF	Price SU
HemCon ChitoGauze Trainer	1 piece	N7 308	24,50



Lifeguard Compact gauze bandages non-sterile

White, z-folded, 6ply, 10 cm x 3,75 m

Vacuum packed, non sterile – for training purposes only

Type	SU	REF	Price SU
non-sterile *	1 piece	J8 1000	1,95

* no medical device



Wound Packing Trainer

This affordable and realistic Wound Packing Simulator is ideal for emergency training and a simulation patient program. Shown is a realistic bullet wound with bleeding function caused by a large gun caliber, ideal for training the Wound Packing procedure. The Simulator is made of a high quality and very durable silicone that simulates human skin and tissue. Common hemostatic wound pads can be used. The trainer is fast and easy to clean after use.

The package contains:

Wound Packing Simulator (Ref.no. 7770)

250 ml artificial blood (Ref.no. 7024)

Bottle for fake blood (Ref.no. 7770-1)

Tube 80 cm (Ref.no. 7770-2)

Plastic tray (Ref.no. 7770-3)

Technical data: Dimensions: 12.5 x 10 x 5 cm, Weight: 1 kg

Type	SU	REF	Price SU
Wound Packing Trainer	1 piece	N7 7770	329,00



Trauma Leg Task Trainer

Practice first aid measures and manage a major hemorrhage with the Trauma Leg. The trainer represents an injured leg with laceration, bullet entry, and bullet exit wounds. You can apply direct pressure, pack the wound with gauze, place a pressure dressing, and apply the tourniquet to really stop the bleeding (tourniquet not included).

- » High Quality Silicone represents soft tissues
- » The bleeding system enables the flow of blood from the wounds
- » Tourniquet application stops the bleeding
- » Simulated blood included
- » Application bottle included

Technical data: Length: 48 cm, Diameter: 18 cm

Type	SU	REF	Price SU
Trauma Leg Task Trainer	1 piece	N7 X1000	1.225,00



SU = Sales unit. All prices in Euro, plus value added tax.



Body Bag ECO

Opaque body bag, waterproof and tearproof, with welded polyester zipper, biodegradable. The cover does not replace a transport device.

Technical data: Colour: white, Material thickness: 0.16 mm, Dimensions: L 2400 x W 900 mm

Type	SU	REF	Price SU
Body Bag Eco	1 piece	N6 1102	13,50



Body bag Standard

Body bag with double round closure in C shape, reinforced with 4 handles, carrying capacity with proper use up to 180 kg (the body bag does not replace a transport device); combustible according to VDI 3891, disinfectant resistant and liquid tight and limited gas-tight.

Technical data: Colour: black, Material thickness: 0.18 mm, Dimensions: L 2340 x B 900 mm

Type	SU	REF	Price SU
Body Bag Standard	1 piece	N6 1104	15,50



Transport body bag

Black transport body bag made of vinyl with double zipper on three sides and 6 handles. Thanks to its particular durability, especially suitable for an emergency or can be used for similar situations.

Technical data: Colour: black, Material thickness: 0.40 mm, Dimensions: L 2150 x W 900 mm

Type	SU	REF	Price SU
Transport Body Bag	1 piece	N6 1105	33,50

of:

sender stamp

Customer number

company

street and house number

postcode	city
----------	------

tel. _____ fax _____e-mail

date and signature

Please enclose the following catalogs with my order:

Your emergency specialist catalog

C

Rapid diagnostic for the practice, laboratory, etc.

C

Tactical medical equipment

We order for the fastest possible delivery, recognising your terms and conditions

[illegible]



CATALOG RANGE

Special catalogs

1



2



3



4



Our special selection:

- | | | | | |
|---|-------------------------------|------------------------------------------------|-----------|---------------|
| 1 | medical products | The complete catalog | 896 pages | Z1 93-P |
| 2 | schnell-diagnostik | Rapid tests for the practice, laboratory, etc. | 112 pages | Z1 04-07-03 |
| 3 | emergency products | Your emergency specialist catalog | 100 pages | Z1 98-PL-E |
| 4 | tactical care products | Tactical medical equipment | 24 pages | Z1 18-03-02-E |

tactical care products

 *Tactical medical equipment*



(91)Z1 18-03-02-E(10)555-7-0001-2501

Your specialist retailer – because no one can advise you better!